

**Table (15): Results of Stability Test of Tadalafil After 24 hrs**

<b>Time and Temp.</b>	<b>Average of AUC of Tadalafil (80 µg/mL)</b>		<b>Assay%</b>	
	Cialis®	Adam®	Cialis®	Adam®
<b>Standard solution</b>				
<b>24 hr at 25°C</b>	980739	1420564	98.93	100.06
<b>24 hrs at 4 °C</b>	1003364	1420449	101.2	100.5

The stability results within stated limit of range 98%-102% in 48 hrs are listed in tables (16,17) and figures (21,22) .

**Table (16): Results of Stability Test of Esomeprazole After 48 hrs**

<b>Time and Temp</b>	<b>Average of AUC of Esomeprazole. (120µg/mL)</b>		<b>Assay%</b>	
	Nexium®	Pumpinox®	Nexium®	Pumpinox®
<b>Standard solution</b>				
<b>48 hrs in 25 °C</b>	1510442	1696718	98.05	99.07
<b>48 hrs in 4 °C</b>	1512311	1687765	98	99.6