

**Table 22: Intra-day precision and accuracy data for LLOQ samples of candesartan based on the standard calibration curve of day three validation.**

Theo. Conc.	AUC Drug	AUC IS	Ratios	Measured Conc.	Precision %	Accuracy %
10 QC Low	3759	49891	0.075	10.888	4.816	108.88
	3583	49509	0.072	10.076		100.76
	3641	49925	0.073	10.230		102.30
	3836	52316	0.073	10.336		103.36
	3716	50837	0.073	10.273		102.73
	3637	52123	0.070	9.368		93.68
Mean	3695.333	50766.833	0.073	10.195		101.952
STD	92.996	1208.465	0.002	0.491		4.910
CV%	2.517	2.380	2.470	4.816		4.816

**Table 23: Intra-day precision and accuracy data for QC low samples of candesartan based on the standard calibration curve of day three validation.**

Theo. Conc.	AUC Drug	AUC IS	Ratios	Measured Conc.	Precision %	Accuracy %
30 ng/ml QC Low	8033	51565	0.156	32.878	5.080	109.59
	7342	50607	0.145	29.949		99.83
	7458	50743	0.147	30.470		101.57
	6694	48096	0.139	28.339		94.46
	6617	46299	0.143	29.356		97.85
	8311	55992	0.148	30.865		102.88
Mean	7409.167	50550.333	0.146	30.310		101.032
STD	685.326	3309.961	0.006	1.540		5.132
CV%	9.250	6.548	3.846	5.080		5.080