

Table 3.28 Inter-day accuracy & precision data for the three days of QC samples validation (LLOQ, QC low, QC mid, and QC high) of propranolol based on the standard calibration curves for the three days of validation (n=6).

	50			150			1500			2500		
	Day One	Day Two	Day Three	Day One	Day Two	Day Three	Day One	Day Two	Day Three	Day One	Day Two	Day Three
	51.764	54.071	50.619	152.764	152.134	156.355	1530.978	1516.459	1454.968	2624.655	2552.107	2556.592
	53.529	53.444	51.528	149.247	156.355	151.082	1552.438	1497.534	1495.227	2592.366	2512.323	2517.108
	52.588	55.075	51.978	147.878	153.268	156.260	1551.392	1491.927	1513.767	2621.938	2532.287	2410.762
	52.488	52.532	50.843	148.955	147.411	154.590	1529.660	1489.235	1527.823	2605.912	2558.965	2591.729
	52.117	52.388	51.146	147.545	150.675	155.216	1544.492	1509.337	1505.963	2579.148	2529.755	2610.097
	50.080	54.865	50.911	149.290	151.834	148.128	1554.274	1524.310	1527.556	2639.846	2543.392	2557.940
Mean	52.331			151.610			1517.630			2563.162		
STD	1.426			3.184			26.016			54.055		
CV%	2.724			2.100			1.714			2.109		
Accuracy %	104.66			101.07			101.18			102.53		