

Table 3.5 Inter-day accuracy & precision data for the three days of QC samples validation (LLOQ, QC low, QC mid, and QC high) of propranolol based on the standard calibration curves of the three days of validation (n=6).

	50			150			1500			2500		
	Day One	Day Two	Day Three	Day One	Day Two	Day Three	Day One	Day Two	Day Three	Day One	Day Two	Day Three
	53.458	53.391	56.198	160.876	161.496	159.757	1619.388	1635.247	1634.039	2611.061	2631.296	2650.168
	54.009	54.632	49.209	164.455	165.704	159.795	1636.410	1475.932	1573.784	2601.889	2640.738	2589.078
	51.689	51.716	54.682	162.866	154.017	167.625	1588.879	1575.259	1605.984	2629.035	2647.646	2581.866
	51.596	52.526	53.511	154.524	156.583	151.162	1626.939	1601.593	1603.539	2592.341	2648.329	2635.382
	51.320	53.355	49.130	148.731	152.624	150.396	1605.938	1584.088	1583.558	2525.086	2484.763	2456.730
	48.075	51.455	52.850	148.385	148.205	150.304	1584.630	1614.234	1620.212	2619.954	2622.595	2568.565
Mean	52.378			156.528			1598.314			2596.473		
STD	2.099			6.427			36.822			56.141		
CV%	4.008			4.106			2.304			2.162		
Accuracy %	104.76			104.35			106.55			103.86		