

HPLC method was developed, validated as per ICH guideline and performed for the determination of esomeprazole magnesium trihydrate in tablets using a C18 column, the mobile phase was ACN/phosphate buffer (60:40,v/v, pH 7) with a flow rate of 1.0 ml/min was applied. UV detection at 205 nm using lansoprazole as an internal standard. The calibration curve of esomeprazole was linear in the range of 100~1000 µg/ml ( $r = 0.9992$ ,  $n=4$ ). And The mean recovery for esomeprazole from tablets ranged between 97.82~98.22% (Onal A., *et al.*, 2006).

As can be seen from the above studies no method was developed for simultaneous determination of esomeprazole and tadalafil in pharmaceutical formulation accordingly we conduct this research.

### **1.8 Literature Survey for Determination of Tadalafil**

RP-HPLC method for estimation of tadalafil in dosage forms was developed and validated , using a Agilent Eclipse XBD C18 column (150 X 4.6mm i.d., 5µm particle size) with mobile phase, acetonitrile and buffer solution (50:50v/v), at a flow rate of 1.2ml/min. at 282nm, methanol was used as diluent (Alivelu Samala,*et al.*,2013).

Direct reverse phase high performance liquid Chromatography method has been published for the estimation of tadalafil in Tablet dosage form . the method was validated and conducted by using a mobile phase composition of (70:30) buffer : acetonitrile, pH adjusted to 4, flow rate and UV detection was 0.8ml/min, 277 nm respectively (Kamepalli Sujana, *et al.*, 2012).

Simultaneous determination of tadalafil and dapoxetine in solid dosage form was achieved through RP-HPLC method, Water Symmetry C-18 (150x4.6mm), 5µ and a mobile phase composed of Buffer : Acetonitrile (65:35). The retention time of Tadalafil